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Accessoires
Dispositif description
Le Teosyal®Pen
The Teosyal®Pen is a new wireless and motorized injection system powered by a lithium battery.
The Teosyal®Pen is a non-sterile device intended to assist medical practitioners in injecting dermal fillers.
After inserting 1ml syringe into the Syringe Holder (TeoLink), the practitioner simply presses the top part of the TeoLink for a smooth and even delivery of the product, in continuous flow mode or in drops.
The box of Teosyal®Pen does not include a battery cap or TeoLinks.
The Teosyal®Pen is only intended for Dermal Filler injections.

Le TeoLink
The TeoLink is a syringe holder that has been designed to improve syringe handling ergonomics, allowing the user to have greater freedom to operate during injections.
The TeoLink's cover touch system provides touch activation along the entire upper surface of the syringe holder.
The TeoLink is compatible with most 1ml Hyaluronic Acid syringes on the market.
The TeoLink is compatible with either needles or cannulas.
The TeoLink is compatible exclusively with Teosyal®Pen.
Do not touch the key board while inserting the battery cap to the Teosyal®Pen.
The TeoLink box contains 25 TeoLink units + 1 battery cap.

Advantages of the Teosyal®Pen
For the patient :
• Significantly less pain
• Highly precise injections
• Reduced side effects

For the practitioner :
• User-friendly & ergonomic design for more comfort
• Easy to use & intuitive system
• Optimized control of the procedure : depth of the injection, homogeneous delivery and precise dosage.

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• Do not dispose of the Teosyal®Pen, but rather, return it to Teoxane SA, stating the reasons for the return.
• Do not use the Teosyal®Pen near electronic equipment which emits electromagnetic fields.

Accessories
Battery cap
Button Battery Lithium 6V /170 mAh
The Teosyal®Pen must be used with a certified battery IEC60086-4 only.

TeoLink
TeoLink (syringe holders) : Plastic PC Makrolon 2458

Cleaning
After each use, the Teosyal®Pen must be cleaned with an adequate cleaning agent before re-using it. Do not immerse the device or its components in water.
We recommend using an appropriate cleaning product in accordance with your local regulation.

Handling conditions
It is recommended to store and transport the Teosyal®Pen in its original packaging. The Teosyal®Pen box can also be used as a display.
Operating conditions 0 - 2000 m
Temperature : -20 °C TO 25 °C
Humidity : 60% RH AT 30 °C non-condensing

Transport conditions
Temperature : -20 °C TO 40 °C
Humidity : 60% RH AT 30 °C non-condensing

Storage condition
Temperature : -20 °C TO 40 °C
Humidity : 60% RH AT 30 °C non-condensing

Assembly and use
Please refer to the assembly instructions drawings attached.
1. Unplug the black button of the Teosyal®Pen
2. Insert the disposable cap without touching the key board of the Teosyal®Pen

3. Screw the battery cap
4. Unscrew the plunger rod from the syringe
5. Discard the plunger rod
6. Remove the finger grip from the syringe
7. Discard the finger grip
8. Insert the syringe into the TeoLink
9. Connect the TeoLink to the Teosyal®Pen
10. "Click" indicates the proper attachment of the TeoLink to the Teosyal®Pen
11. Remove and discard the cap from the syringe
12. Attach the needle to the syringe
13. Remove the needle protection. Do not discard it
14. Press the ON/OFF button. The green light indicates that the device is switched ON
15. Select the delivery mode and set-up the delivery speed or drop size
16. By pressing the top of the TeoLink move the lead screw moves forward until it reaches the plunger
17. The injection can be performed by pressing the top of the TeoLink

After the use
18. Replace the needle protection
By pressing on two sides of the TeoLink release the syringe holder with the syringe and needle
20. Discard the TeoLink, the syringe and the needle
21. Press the reverse button to return the lead screw to its original position

The Device Control Buttons and Control Lights
1. Activation Button
2. Device ON/OFF button and light indicator
3. Continuous flow mode speed set-up button and light indicator
4. Drop mode drop size set-up button and light indicator
5. Original position screw button and light indicator

Warranty
The warranty period for the Teosyal®Pen is 18 months from the date of purchase.

Note
Teoxane SA is not responsible for the cost of replacement of unused injectable material or injectable material lost or wasted during set-up, operation or maintenance of the system.
Do not use or contact Teoxane SA or your local representative for any further questions.
• The syringe holders (TeoLink), which are designed for use with the Teosyal®Pen exclusively and therefore cannot be used with any other third party device.
• Do not use or dispose of any part of the Teosyal®Pen or TeoLink. Damage of the Teosyal®Pen caused by third party batteries voids the warranty.
• Do not sterilize any of the components : Teosyal®Pen, TeoLink (syringe holder) and battery.
• Wearing single-use gloves is highly recommended for use of the equipment.
• No alteration of this equipment is allowed.
• Do not use the TeoLink (syringe holder) or battery according to your local environmental regulations.

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